



PRODUCT CODE: 272170

Buffer Solution pH 7.00 (20°C)

TARIC 3822 00 00 00

The measure of pH in the quality control of our buffer solutions is traceable with the formulation according to U.S. National Institute of Standards and Technology Buffers .

PHYSICAL DATA: liquid, Miscible with water • D 20/4 1,01 • pH7,0 •

BIBLIOGRAPHY:

HAZARDOUS:

WEIGHT/VOLUME INFORMATION: 1l~1,006 kg 1kg~0,994 l

COMPOSITION:

Potassium di-Hydrogen Phosphate	3,522 g
di-Sodium Hydrogen Phosphate 12-hydrate	14,020 g
Preservative 1%	5 ml
Water (s.q.m)	1 l

SPECIFICATIONS:

pH VALUES DEPENDING ON TEMPERATURE:	
T (°C)	ApH

0	+0,13
5	+0,07
10	+0,05
15	+0,02
20	
25	-0,02
30	-0,02
35	-0,04
40	-0,05
50	-0,05
Uncertainty	See certificate
Traceability	NIST
pH at 20°C	6,98 - 7,02

Ed.: 2 . Vig.: 10.07.2014 .

Prod.: 272170