

# TECHNICAL DATA SHEET

---

<b>Part Number</b>	55674-50G-HON
<b>Description</b>	Ammonium formate

---

## Technical specifications :

- ✓ ammonium formate, eluent additive for lc-ms 50g
- ✓ CAS N°540-69-2
- ✓ Chemical purity ≥99%
- ✓ Molecular weight 63.06 g/mol
- ✓ Molecular formula CH<sub>5</sub>NO<sub>2</sub>
- ✓ Density 1.266 g/cm<sup>3</sup> (20 °C)
- ✓ Flash point Not applicable
- ✓ Material form crystalline
  
- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Bottle
- ✓ Package size 50g
  
- ✓ GHS pictograms GHS07
- ✓ H Statements H319
- ✓ P Statements P280

## Sales unit :

50g Ammonium formate

---