

Standard Analysis Certificate

PRODUCT :CONDUCTIVITY STANDARD 200 µS/cm RS
CODE :575070
METHOD :10454

| TEST | U.M. | SPECIFICATION |
|----------------------|-------|-------------------------|
| Description | - | Clear colourless liquid |
| Identification | - | Positive |
| Conductivity at 25°C | µS/cm | 198.0 ÷ 202.0 |

Traceable to Standard NIST N° 999a
 Confidence interval: ± 1 µS/cm

DOMINIQUE DUTSCHER SAS

Date :11/01/2020

QUALITY CONTROL RESPONSIBLE
B. COULANGE (VDR)