

Standard Analysis Certificate

PRODUCT :CONDUCTIVITY STANDARD 500 µS/cmRS
CODE :575080
METHOD :10455

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Identification	-	Positive
Conductivity at 25°C	µS/cm	495.0 ÷ 505.0

Traceable to Standard NIST N° 999a
 Confidence interval: ± 1 µS/cm

DOMINIQUE DUTSCHER SAS

Date :11/01/2020

QUALITY CONTROL RESPONSIBLE
B. COULANGE (VDR)