

Standard Analysis Certificate

PRODUCT :CONDUCTIVITY STANDARD 20000 μ S/cm RS
CODE :575140
METHOD :10461

TEST	U.M.	SPECIFICATION
Description	-	Clear colourless liquid
Identification	-	Positive
Conductivity at 25°C	μ S/cm	19800 \div 20200

Traceable to Standard NIST N° 999a
 Confidence interval \pm 100 μ S/cm

DOMINIQUE DUTSCHER SAS

Date :11/01/2020

QUALITY CONTROL RESPONSIBLE
B. COULANGE (VDR)