


PDS No. 656197	PRODUCT DATA SHEET	Page 1 of 1
Revision 01	Lid for Microplates, PS	 greiner bio-one
	Greiner Item-No. 656197	

1.	Description / Specification	
1.1	Description	PS lid for Microplate, single position, high profile, without condensation rings
1.2	Dimensions	Length: 127,5 mm (+/- 0,2 mm) Width: 85 mm (+/- 0,2 mm) Curvature: +/- 300 µm
1.3	Volume	N/A
1.4	Material / Resin	PS (Polystyrene), free of heavy metal
1.5	Colour	Clear
1.6	Sterilization	No
1.7	Quality Control	- Raw Material-Control: physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	-
2.2	Temperature range	-20°C to +60°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: www.gbo.com/bioscience →Products →Literature →Technical Information →Chemical Resistance of Resins
2.6	Shelf life	N/A
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	10
3.2	Pieces / Box	150
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	-

4.	Other Information	

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision -	Date 17 December 2010	Date 17 December 2010	Date 17 December 2010	
Date -	Name S. Kaelberer	Name Dr. R. Heller	Name A. Schulz	