


PDS No. 662145	PRODUCT DATA SHEET	Page 1 of 1
Revision 01	CrystalBridge™	 greiner bio-one
	Greiner Item-No. 662 145	

1.	Description / Specification	
1.1	Description	CrystalBridge™ insert, concave well bottom for sitting drop experiments, fitting exactly into the wells of the ComboPlates™ (Cat. No. 662050 / 662150 / 662850)
1.2	Dimensions	∅ Well: 4,6 mm
1.3	Volume	Max. volume: 45 µl
1.4	Material / Resin	PS (Polystyrene), free of heavy metal
1.5	Colour	Clear
1.6	Sterilization	No
1.7	Quality Control	<u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	-
2.2	Temperature range	-20°C to +60°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: www.gbo.com/bioscience →Products →Literature →Technical Information→Chemical Resistance of Resins
2.6	Shelf life	N/A
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	250
3.2	Pieces / Box	Quantity on request
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	-

4.	Other Information	
		-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision -	Date 4 February 2010	Date 2 December 2010	Date 3 December 2010	
Date -	Name S. Kaelberer	Name Dr. U. Honisch	Name A. Schulz	