


PDS No. 662825	PRODUCT DATA SHEET		Page 1 of 1
Revision 01	24 Well Intermediate Lid		
	Greiner Item-No. 662825		
Valid for Item-No.:	662825 (sterile)		

1.	Description / Specification	
1.1	Description	24 Well Intermediate Lid, Sterile, for magnetic levitation with 24 Well Bio-Assembler Kit
1.2	Dimensions	See Customer Drawing
1.3	Volume	N/A
1.4	Material / Resin	PET (Polyethylene Terephthalate), free of heavy metal
1.5	Colour	Clear
1.6	Sterilization	SAL 10 ⁻³
1.7	Quality Control	<u>Product-Control:</u> testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for covering or capping to be used by qualified personnel in a laboratory environment.
1.9	Other Information	For single use only

2.	Features	
2.1	Basic features	Contents non-cytotoxic
2.2	Temperature range	N/A
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	5 years
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	2
3.2	Pieces / Box	10
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	-

4.	Other Information

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 00	Date 3 February 2023	Date 12 April 2023	Date 13 April 2023	
Date 06.07.2020	Name S. Kaelberer	Name P. Wachter	Name Dr. Chen-Ket Chai	

DISCLAIMER: The description of a certain product can only be considered as a guidance because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.