

TECHNICAL DATA SHEET

Part Number	60221-250G-HON
Description	Potassium phosphate monobasic

Technical specifications :

- ✓ potassium phosphate monobasic, for hplc =99,5% 250g
- ✓ CAS N°7778-77-0
- ✓ Chemical purity ≥99.5%
- ✓ Molecular weight 136.09 g/mol
- ✓ Molecular formula H₂KO₄P
- ✓ Density 2.34 g/cm³ (20 °C)
- ✓ Flash point Not applicable
- ✓ Material form crystalline

- ✓ Storage condition Room Temperature
- ✓ Package type Glass Bottle
- ✓ Package size 250g

- ✓ P Statements P280

Sales unit :

250g Potassium phosphate monobasic
