

# HJ-BIOANALYTIK

Mikrotiterplatten - Technologie

## Certificate of Quality

Part-no.: 606200  
Lot-no.: 230080  
Description: 96 Micro Tube Rack System  
Type: with tube strips with 8 tubes  
Volume: 0,65 ml  
Material: Polypropylene  
Size: SBS-Standard  
Quantity: 10 racks per pack  
Range of temperature: -80°C up to +100°C  
Chemical resistance: High resistance against biological and organic solvents

This product is free of DNase, RNase and pyrogene.

Shelf-Life: Minimum 10 years  
Origin: Germany  
Manufacturer: HJ-BIOANALYTIK GmbH  
Ferdinand-Clasen-Str. 34  
41812 Erkelenz  
Phone: +49 2431 948415  
Fax: +49 2431 948416  
email: [info@hj-bioanalytik.de](mailto:info@hj-bioanalytik.de)

QC – Manager: Sascha Kuffner



D-Erkelenz, 15.01.2025

.....  
Signature