


## Product specification

**Product number : 61.100.24**

	Specifications
Client	Surpodo
Dimensions slides	75.8 – 75.3mm and 25.8 – 25.5mm
Lay-out slide	 <p><b>61.100.24</b>                      1 well 18 mm</p>
Color	Black
Description	1 well of 18 mm
Material slide :	Soda Glass with frosted end (20 mm)
Labelling	1 label per box, mentioning : <ul style="list-style-type: none"><li>- Product number</li><li>- Product description + color</li><li>- LOT number + number of slides in a box</li><li>- Manufactured by : Immuno-Cell Int. Bvba</li></ul>
Packaging	100 slides/box 10 boxes in outer white box
Quality	Free of foreign matter, dirt, oil and grease + free of visible marks and scratches
Storage and shipment conditions	Ambient temperature : above 7 °C

Rev.0 on 20/01/2020 – Approved by Hanne Vermeeren