	61.100.71
---	------------------

Client	Immuno-Cell
Description	18 wells - 6 mm diameter

PRODUCTION	
Glass size	76 x 26 mm
Frosted end	15 mm
Print	Black
Surface	free of foreign matter, dirt, oil and grease + free if visible marks and scratches


LABELING and PACKAGING	
Carton	100 slides in a box
Labeling	Lotnumber + reference number Immuno-Cell + name + # slides per box

CERTIFICATE	NA
--------------------	----

SHIPMENT CONDITIONS	Ambient Temperature
----------------------------	---------------------

STORAGE CONDITIONS	Ambient Temperature
---------------------------	---------------------

DOCUMENTS REQUIREMENTS	NA
-------------------------------	----



61.100.71

18 wells 6 mm

DOMINICUS E. DUTSCHER STAS