


Product specification

Product number : 61.100.75

	Specifications
Client	Immuno-Cell
Dimensions slides	75.8 – 75.3mm and 25.8 – 25.5mm
Lay-out slide	 61.100.75 4 wells 7 mm
Color	Black
Description	4 wells – 7 mm - numbered
Material slide :	Soda Glass with frosted end (20 mm)
Material number :	62.243.00
Screen Type :	77S
Labelling	1 label per box, mentioning : <ul style="list-style-type: none">- Product number- Product description + color- LOT number + number of slides in a box- Manufactured by : Immuno-Cell Int. Bvba
Packaging	100 slides/box 10 boxes in outer white box
Quality	Free of foreign matter, dirt, oil and grease + free of visible marks and scratches
Storage and shipment conditions	Ambient temperature : above 7 °C

Rev.1 on 15/11/2019 – Approved by Hanne Vermeeren

Immuno-Cell International BVBA - Motstraat 72/8 - B-2800 Mechelen - Belgium
Info@immuno-cell.com - www.immuno-cell.com

Page 1 of 1