


PDS No. 614101	PRODUCT DATA SHEET	Page 1 of 1
Revision 03	Macro Cuvette, 4 ml	 greiner bio-one
	Item-No. 614101	

1.	Description / Specification	
1.1	Description	Macro Cuvette, 4 ml
1.2	Dimensions	See Customer Drawing
1.3	Volume	Working volume: 2,5 ml
1.4	Material / Resin	PS (Polystyrene), free from heavy metal
1.5	Colour	Clear
1.6	Sterilization	No
1.7	Quality Control	<u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	Applicable wavelength range: 340 nm to 900 nm
2.2	Temperature range	-20°C to +60°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	N/A
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Inner Box	100 / Styrofoam package
3.2	Pieces / Box	1000
3.3	Lot-No.	-
3.4	Other Information	-

4.	Other Information	
		-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 02	Date 17 April 2015	Date 20 April 2015	Date 20 April 2015	
Date 13.03.2013	Name S. Kaelberer	Name Dr. T. Schreiber	Name A. Schulz	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.