


PDS No. 6278xx	<b>PRODUCT DATA SHEET</b>			Page 1 of 1
Revision 05	CELLview™, Cell Culture Dish, 35 x 10 mm, Glass Bottom, Sterile			
	Greiner Item-No. 6278xx			
Valid for Item-No.:	<b>627860</b> (sterile)	<b>627861</b> (sterile)	<b>627870</b> (sterile)	<b>627871</b> (sterile)

1.	Description / Specification	
1.1	Description	Cell culture dish with glass bottom for microscopic applications, with vents, sterile <u>Glass bottom</u> : maximal spectral transmission, no autofluorescence, no depolarisation of light 627860: with one compartment, physical surface treatment 627861: with one compartment, without treatment 627870: with four compartments, physical surface treatment 627871: with four compartments, without treatment
1.2	Dimensions	See Customer Drawing <u>Weight</u> : 627860, -861: lid: 1,2 – 1,5 g / dish: 1,6 – 1,8 g 627870, -871: lid: 1,2 – 1,5 g / dish: 2,5 – 2,8 g <u>Glass thickness</u> : 0,17 mm (+/- 0,01 mm)
1.3	Volume	627860, -861: working volume: 5,0 ml 627870, -871: working volume: 0,1-0,5 ml / per compartment 627860, -861: growth area: 8,7 cm <sup>2</sup> 627870, -871: growth area: 1,9 cm <sup>2</sup> / per compartment
1.4	Material / Resin	<u>Dish</u> : PS (Polystyrene), free of heavy metal <u>Bottom</u> : borosilicate glass; hydrolytic class 1 (DIN-ISO 719) <u>Lid</u> : PS (Polystyrene), free of heavy metal
1.5	Colour	<u>Dish</u> : clear <u>Glass Bottom</u> : high transparent <u>Lid</u> : clear
1.6	Sterilisation	SAL 10 <sup>-3</sup>
1.7	Quality Control	- <u>Raw Material-Control</u> : physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	- For single use only - Manufactured according to ISO 8255-1:2011

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic <u>Adhesive</u> : resistant to DMSO and 70% ethanol, compatible with cell culture media
2.2	Temperature range	+ 4°C to +37°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: <a href="https://www.gbo.com/en_INT/know-how-services/download-center.html">https://www.gbo.com/en_INT/know-how-services/download-center.html</a>
2.6	Shelf life	2 years after month of production
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	10
3.2	Pieces / Box	40
3.3	Lot-No.	E JJ MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	Certificate of Quality

4.	Other Information
4.1	Research use only. Not for diagnostics.

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 04	Date 7 July 2015	Date 8 July 2015	Date 8 July 2015	
Date 03.03.2015	Name S. Kaelberer	Name Dr. L. Breth	Name Dr. A. Ganser	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.