


PDS No. 6289x0	<b>PRODUCT DATA SHEET</b>				Page 1 of 1
Revision 02	Cell Culture Dish, CELLCOAT <sup>®</sup> , 60 x 15 mm				 <b>greiner bio-one</b>
	Greiner Item-No. 6289x0				
Valid for Item-No.:	<b>628910</b>	<b>628920</b>	<b>628930</b>	<b>628940</b>	<b>628950</b>

1.	Description / Specification	
1.1	Description	Cell Culture Dish, CELLCOAT <sup>®</sup> , 60 x 15 mm, with vents 628910: Laminin protein coating 628920: Fibronectin protein coating 628930: Poly-L-Lysine protein coating 628940: Poly-D-Lysine protein coating 628950: Collagen Type I protein coating
1.2	Dimensions	Dish: Ø (outside): 54,39 mm (+0,1/-0,2 mm) height: 13,8 mm (±0,1mm) Dish with lid: Ø (outside): 59 mm (+0,1/-0,2 mm) height: 15,4 mm (0,2/-0,1 mm) Total weight: 7,4 – 7,8 g
1.3	Volume	Max. volume: 28 ml Working volume: 6-7 ml Growth area: 21 cm <sup>2</sup>
1.4	Material / Resin	Dish: 628910: PS (Polystyrene), coated with Laminin 628920: PS (Polystyrene), coated with Fibronectin 628930: PS (Polystyrene), coated with Poly-L-Lysine 628940: PS (Polystyrene) coated with Poly-D-Lysine 628950: PS (Polystyrene), coated with Collagen Type I Dish / Lid: the materials for manufacturing are free of heavy metals
1.5	Colour	Dish: clear Lid: clear
1.6	Sterilisation	No
1.7	Quality Control	- Raw Material-Control: physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	- For single use only - Elevations for ventilation of culture

2.	Features	
2.1	Basic features	-
2.2	Temperature range	628910, -920: +2°C to +8°C 628930, -940, -950: room temperature
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: <a href="http://www.gbo.com/bioscience">www.gbo.com/bioscience</a> →Products →Literature →Technical Information→Chemical Resistance of Resins
2.6	Shelf life	628910, -920: 6 month after month of production 628930, -940, -950: 2 years after month of production
2.7	Other Information	-

3.	Packaging	628910, -920	628930, -940, -950
3.1	Pieces / Bag	5	20
3.2	Pieces / Box	20	100
3.3	Lot-No.	AABB0D89 (year, week, country code, index, unique code)	AABB0D89 (year, week, country code, index, unique code)
3.4	Other Information	-	-

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 01	Date 15 February 2013	Date 18 February 2013	Date 18 February 2013	
Date 16.03.2011	Name S. Kaelberer	Name Dr. T. Schreiber	Name A. Schulz	