


PDS No. 633181	PRODUCT DATA SHEET	Page 1 of 1
Revision 02	Petri Dish With Vents, Sterile	 greiner bio-one
	Greiner Item-No. 633181	

1.	Description / Specification	
1.1	Description	Petri dish with vents, 94 x 16 mm, sterile
1.2	Dimensions	<u>See Customer Drawing</u> Total weight: 13,0 – 14,7
1.3	Volume	Max. volume: 80 ml Working volume: 15 - 16 ml
1.4	Material / Resin	<u>Dish</u> : PS (Polystyrene), free of heavy metal <u>Lid</u> : PS (Polystyrene), free of heavy metal
1.5	Colour	<u>Dish</u> : clear <u>Lid</u> : clear
1.6	Sterilisation	SAL 10 ⁻³
1.7	Quality Control	- <u>Raw Material-Control</u> : physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	- For single use only - Elevations for ventilation of culture

2.	Features	
2.1	Basic features	-
2.2	Temperature range	-20°C to +60°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	5 years after month of production
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	20
3.2	Pieces / Box	480
3.3	Lot-No.	F YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	-

4.	Other Information	
		-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 01	Date 2 March 2015	Date 3 March 2015	Date 3 March 2015	
Date 24.03.2011	Name S. Kaelberer	Name Dr. T. Schreiber	Name A. Schulz	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.