


PDS No. 6391xx/ 639160-TRI	<b>PRODUCT DATA SHEET</b>			Page 1 of 1
Revision 11	Petri Dish, 145 x 20 mm			
	Greiner Item-No. 6391xx / 639160-TRI			
Valid for Item-No.:	<b>639102</b>	<b>639160</b> (sterile)	<b>639161</b> (sterile)	<b>639160-TRI</b> (sterile)

1.	Description / Specification	
1.1	Description	Petri Dish with vents 639102: standard 639160: CELLSTAR <sup>®</sup> , physical surface treatment, sterile 639161: sterile 639160-TRI: CELLSTAR <sup>®</sup> , physical surface treatment, sterile, triple packed
1.2	Dimensions	See Customer Drawings Weight: 639102, -161: 42,0 – 46,0 g (total weight) 639160, -160-TRI: 21,1-22,5 g (lid) 24,0-25,4 g (dish)
1.3	Volume	Max. volume: 240 ml Working volume: 25 - 27 ml 639160, -160-TRI: growth area: 143 cm <sup>2</sup>
1.4	Material / Resin	Dish: PS (Polystyrene), free of heavy metal Lid: PS (Polystyrene), free of heavy metal
1.5	Colour	Dish: clear Lid: clear
1.6	Sterilisation	639102: no 639160, -161: SAL 10 <sup>-3</sup> 639160-TRI: SAL 10 <sup>-6</sup>
1.7	Quality Control	Raw Material-Control: physical testing Product-Control: testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	- For single use only - Elevations for ventilation of culture

2.	Features	
2.1	Basic features	639102, -161: - 639160, -160-TRI: free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2	Temperature range	-20°C to +60°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: <a href="https://www.gbo.com/en_INT/know-how-services/download-center.html">https://www.gbo.com/en_INT/know-how-services/download-center.html</a>
2.6	Shelf life	639102, -161: 5 years after month of production 639160, -160-TRI: 4 years after month of production
2.7	Other Information	-

3.	Packaging	639160	639102, -161	639160-TRI
3.1	Pieces / Bag	5	15	5 (triple packed)
3.2	Pieces / Box	120	120	120
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	F YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	Certificate of Quality	-	Certificate of Quality

4.	Other Information
4.1	639160 / 639160-TRI: Research use only. Not for diagnostics.

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 10	Date 31 July 2015	Date 3 August 2015	Date 26 August 2015	
Date 03.03.2015	Name S. Kaelberer	Name Dr. L. Breth	Name Dr. A. Ganser	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.