


PDS No. 653081	<b>PRODUCT DATA SHEET</b>	Page 1 of 1
Revision 02	Microscopic Cover Glass	 <b>greiner bio-one</b>
	Greiner Item-No. 653081	

1.	Description / Specification	
1.1	Description	Microscopic Cover glass is used in medical and biological research as cover material for specimens in microscopic examinations. Suited for microplates.
1.2	Dimensions	Length: 76 mm Width: 51 mm Thickness: 0,145 mm (+/-0,015 mm)
1.3	Volume	-
1.4	Material / Resin	Borosilicate glass
1.5	Colour	Clear
1.6	Sterilization	No
1.7	Quality Control	<u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	High chemical resistance
2.2	Temperature range	N/A
2.3	Autoclavability	Yes
2.4	Centrifugation, max. RCF	-
2.5	Chemical Resistance	Hydrolytic resistance acc. to DIN ISO 719 Acid resistance acc. to DIN 12116 Alkali resistance acc. to DIN ISO 695
2.6	Shelf life	N/A
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	50
3.2	Pieces / Box	50
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	-

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 01	Date 26 November 2014	Date 28 November 2014	Date 28 November 2014	
Date 01.06.2012	Name S. Kaelberer	Name Dr. U. Honisch	Name A. Schulz	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.