


PDS No. 6581xx / 6581xx-TRI	PRODUCT DATA SHEET			Page 1 of 1
Revision 10	Cell Culture Flask, CELLSTAR [®] , 250 ml			
	Greiner Item-No. 6581xx / 6581xx-TRI			
Valid for Item-No.:	658170 (sterile)	658175 (sterile)	658190 (sterile)	658195 (sterile)
	658170-TRI (sterile)	658175-TRI (sterile)	658190-TRI (sterile)	658195-TRI (sterile)

1.	Description / Specification	
1.1	Description	Cell Culture Flask, 250 ml, canted neck, printed graduation and labelling field on the side (scale 20-125 ml), sterile 658170: physical surface treatment, standard screw cap* 658170-TRI: physical surface treatment, standard screw cap*, triple packed 658175: physical surface treatment, filter cap 658175-TRI: physical surface treatment, filter cap, triple packed 658190: hydrophobic surface for suspension culture, standard screw cap* 658190-TRI: hydrophobic surface for suspension culture, standard screw cap*, triple packed 658195: hydrophobic surface for suspension culture, filter cap 658195-TRI: hydrophobic surface for suspension culture, filter cap, triple packed * with ventilation position
1.2	Dimensions	Flask: see Customer Drawing 658175, 175-TRI, -195, -195-TRI: pore size of filter membrane: 0.2 µm
1.3	Volume	Max. volume: 250 ml 658170, -170-TRI, -175, -175-TRI: - working volume: 15-38 ml - growth area: 75 cm ²
1.4	Material / Resin	Flask: PS (Polystyrene), free of heavy metal Cap: HDPE (High Density Polyethylene), free of heavy metal Filter: PET (Polyethylene Terephthalate), free of heavy metal
1.5	Colour	Flask: clear; Print: white 658170, -170-TRI, -175, -175-TRI: Cap: red 658190, -190-TRI, -195, -195-TRI: Cap: white 658175, -175-TRI, -195, -195-TRI: Filter: white
1.6	Sterilization	658170, -175, -190, -195: SAL 10 ⁻³ 658170-TRI, -175-TRI, -190-TRI, -195-TRI: SAL 10 ⁻⁶
1.7	Quality Control	- Raw Material-Control: physical testing - Product-Control: testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	- For single use only - Expiry date and Lot-No. printed on bottom of flask 658190, -190-TRI, -195, -195-TRI: additional "Suspension" printed on bottom of flask

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2	Temperature range	-20°C to +60 °C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	658170, -170-TRI, -175, -175-TRI: 4 years after month of production 658190, -190-TRI, -195, -195-TRI: 5 years after month of production
2.7	Other Information	-

3.	Packaging	658170, -175, -190, -195	658170-TRI, -175-TRI, -190-TRI, -195-TRI
3.1	Pieces / Bag	5	5 (triple packed)
3.2	Pieces / Box	120	90
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	
3.4	Other Information	Certificate of Quality	
4.	Other Information		
4.1	Laboratory use only		

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 09	Date 30 November 2016	Date 20 June 2017	Date 21 June 2017	
Date 18.08.2016	Name S. Kaelberer	Name Dr. L. Breth	Name A. Schulz	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.