


PDS No. 6589x0	<b>PRODUCT DATA SHEET</b>				Page 1 of 1
Revision 06	Cell Culture Flask, CELLCOAT <sup>®</sup> , 75 cm <sup>2</sup>				 <b>greiner bio-one</b>
	Greiner Item-No. 6589x0				
Valid for Item-No.:	<b>658910</b>	<b>658920</b>	<b>658940</b>	<b>658950</b>	

1.	Description / Specification	
1.1	Description	Cell Culture Flask, CELLCOAT <sup>®</sup> , 75 cm <sup>2</sup> , canted neck, printed graduation (scale 20-125 ml), writing area with filter screw cap 658910: Laminin protein coating 658920: Fibronectin protein coating 658940: Poly-D-Lysine protein coating 658950: Collagen Type I protein coating
1.2	Dimensions	Flask: see Customer Drawing Pore size of filter membrane: 0,2 µm
1.3	Volume	Max. volume: 250 ml Working volume: 15-38 ml Growth area: 75 cm <sup>2</sup>
1.4	Material / Resin	Flask: 658910: PS (Polystyrene) coated with Laminin 658920: PS (Polystyrene) coated with Fibronectin 658940: PS (Polystyrene) coated with Poly-D-Lysine 658950: PS (Polystyrene) coated with Collagen Type I Screw cap: HDPE (High Density Polyethylene) Filter: PET (Polyethylene Terephthalat), PTFE (Polytetrafluorethane) The materials for manufacturing are free of heavy metals
1.5	Colour	Flask: clear; print: graduation and writing area white Screw cap: red; filter: white
1.6	Sterilization	No
1.7	Quality Control	- <u>Raw Material-Control</u> : physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	-
2.2	Temperature range	658910, -920: +2°C to +8°C 658940, -950: room temperature
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: <a href="https://www.gbo.com/en_INT/know-how-services/download-center.html">https://www.gbo.com/en_INT/know-how-services/download-center.html</a> (Only concerning the standard flask without coating)
2.6	Shelf life	658910, -920: 6 month after month of production 658940, -950: 2 years after month of production
2.7	Other Information	-

3.	Packaging	658910, -920	658940, -950
3.1	Pieces / Bag	10	5
3.2	Pieces / Box	10	50
3.3	Lot-No.	AABB0D89 (year, week, country code, index, unique code)	AABB0D89 (year, week, country code, index, unique code)
3.4	Other Information	-	-

4.	Other Information
4.1	Research use only. Not for diagnostics.

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 05	Date 20 April 2015	Date 21 April 2015	Date 21 April 2015	
Date 02.12.2014	Name S. Kaelberer	Name Dr. T. Schreiber	Name A. Schulz	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.