


| | | |
|---|---|------------------------------------|
|  | ANNEX (AN) | G AN QMI 02.01.01-01 Rev.22 |
| | Medical Products and In Vitro Diagnostic Devices | Seite 1/1 Gültig ab: 21.11.2016 |

¹⁾ in accordance with EU directive 93/42/EC

²⁾ in accordance with EU directive 98/79/EC

| Item-No. | Description | Classification |
|--------------------|--|--|
| | Cryo.s™: | IVD²⁾, CE |
| 121 2xx-xxx | Cryo.s™ 2ml, round bottom, w/o Barcode | |
| 122 2xx-xxx | Cryo.s™ 2ml, round bottom, starfoot base, w/o Barcode | |
| 123 2xx-xxx | Cryo.s™ 1ml, conical bottom, w/o Barcode | |
| 124 2xx-xxx | Cryo.s™ 5ml, round bottom, w/o Barcode | |
| 126 2xx-xxx | Cryo.s™ 2ml, round bottom, starfoot base, w/o Barcode | |
| 127 2xx-xxx | Cryo.s™ 4ml, round bottom, starfoot base, w/o Barcode | |
| | Faeces Container: | IVD²⁾, CE |
| 430 102 | Faeces Container (18.5/61 mm) | |
| 431 102 | Faeces Container (22/63 mm) | |
| 433 102 | Faeces Container (22/63 mm) | |
| 434 102 | Faeces Container (41/57 mm) | |
| 443 102 | Faeces Container (24.9/90 mm) | |
| | Sputum Container: | IVD²⁾, CE |
| 432 102 | Sputum Container (22/63 mm) | |
| | ParoCheck® Collection Kit: | MP Class I¹⁾, CE |
| 460 070 | ParoCheck® Collection Kit | |
| | ParoCheck®: | IVD²⁾, CE |
| 460 010 | ParoCheck® 10, for 10 samples | |
| 460 110 | ParoCheck® 10, for 60 samples | |
| 460 001 | ParoCheck® 10, for 1200 samples | |
| 460 020 | ParoCheck® 20, for 10 samples | |
| 460 120 | ParoCheck® 20, for 60 samples | |
| 460 002 | ParoCheck® 20, for 1200 samples | |
| | PapilloCheck®: | IVD²⁾, CE |
| 465 060 | PapilloCheck®, for 60 samples | |
| 465 088 | PapilloCheck®, for 288 samples | |
| | PapilloCheck® high-risk: | IVD²⁾, CE |
| 505 060 | PapilloCheck® high-risk, for 60 samples | |
| | PapilloCheck® Collection Kit: | |
| 465 070 | PapilloCheck® Collection Kit | MP Class I¹⁾, CE |
| 465 075 | PapilloCheck® Collection Kit | MP Class I¹⁾, CE |
| 862 080 | CheckReport™ Software | Active IVD-accessory²⁾ |
| 862 070 | CheckScanner™ | Active IVD-accessory²⁾ |
| 515 040 | oCheck® DNA Extraction Kit | IVD-accessory²⁾ |
| 517 070 | oCheck® DNA Extraction Kit - CheckExtractor | IVD-accessory²⁾ |
| 517075 | Extension Kit for oCheck® DNA Extraction Kit - CheckExtractor | IVD-accessory²⁾ |
| 865 070 | CheckExtractor™ | IVD-accessory²⁾ |