


Specification

**Protein Marker VI (10 - 245) prestained**

**A8889**

<b>Product Code:</b>	A8889
<b>Product Name:</b>	Protein Marker VI (10 - 245) prestained
<b>Specifications:</b>	<ul style="list-style-type: none"> <li>: • Number of bands: 12</li> <li>: • Three-color protein standard with 12 pre-stained. Proteins are covalently coupled with a blue chromophore except for two reference bands (one green and one red band at 25 kDa and 75 kDa respectively).</li> <li>: • Fragment size range in kDa: 10 - 245</li> </ul>
<b>WGK:</b>	1
<b>Storage:</b>	-20°C
<b>Shipment:</b>	wet ice
<b>CS:</b>	38221900
	

**AppliChem GmbH**

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • [info.de@itwreagents.com](mailto:info.de@itwreagents.com) • [www.itwreagents.com](http://www.itwreagents.com)  
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 07 228 16476 • Register court Darmstadt HRB Nr. 7340