

Specification

Tris for buffer solutions

A1379

Solubility:	800 g/L (H ₂ O)
Physical Description:	Solid
Product Code:	A1379
Product Name:	Tris for buffer solutions
Specifications:	Assay (titr.): min. 99.3 % pH (1 M; H ₂ O; 20°C): 10.5 - 12.0 Heavy metals (as Pb): max. 0.001 % Water: max. 0.5 %
WGK:	1
Storage:	RT
Molecular Formula:	C ₄ H ₁₁ NO ₃
M:	121.14 g/mol
CAS:	77-86-1
EINECS:	201-064-4
CS:	29221900

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode
 DRESDEFF508 • Finanzamt Darmstadt 07 228 16476 • Register court Darmstadt HRB Nr. 7340