

Specification

Water PCR tested, DNA free, for molecular biology, suited for qPCR

A8510

Physical Description:	Liquid
Product Code:	A8510
Product Name:	Water PCR tested, DNA free, for molecular biology, suited for qPCR
Specifications:	DNases/RNases: not detectable Bacterial DNA: not detectable (min. 40 PCR cycles)
WGK:	nwg
Storage:	RT
Molecular Formula:	H ₂ O
M:	18.01 g/mol
CAS:	7732-18-5
EINECS:	231-791-2
CS:	28539010

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 07 228 16476 • Register court Darmstadt HRB Nr. 7340