



**Model:**  
Camera adapter

**Description:**  
CCD camera adapter 0.5x

<b>Optical magnification</b>	0.5x
<b>Connection to the camera</b>	C-mount standard
<b>CS-mount</b>	Allowed, but the C-mount ring (optional) is suggested
<b>Connection to the microscope</b>	23 mm eyepiece connection
<b>Allowed connections</b>	30 mm & 30.5mm eyepiece connection ( with M-113.1 & M-113.2 )
<b>Sensor size compatibility</b>	1/2" 1/2.3" 1/2.5"
<b>External fast identification</b>	Black
<b>Length</b>	55 mm
<b>Color</b>	Black
<b>Allowed</b>	The camera is free to move 360 degree