

Specification

Agarose low EEO (Agarose Standard)

A2114

Physical Description:	Solid
Product Code:	A2114
Product Name:	Agarose low EEO (Agarose Standard)
Specifications:	DNases/RNases/Proteases: not detectable Ash: $\leq 0.5 \%$ Electroendoosmosis (EEO): 0.09 - 0.13 Moisture: $\leq 10 \%$ Gel strength (1 %): $\geq 1200 \text{ g/cm}^2$ Gel strength (1.5 %): $\geq 2300 \text{ g/cm}^2$ Gel point (1.5 %): $36 \pm 1.5^\circ\text{C}$ Melting point: $88 \pm 1.5^\circ\text{C}$ Sulfate: $\leq 0.2 \%$
WGK:	1
Storage:	RT
CAS:	9012-36-6
EINECS:	232-731-8
CS:	39139000

AppliChem GmbH