

Specification

**Dextran Sulfate 500 Sodium Salt *BioChemica***

**A2250**

<b>Physical Description:</b>	Solid
<b>Product Code:</b>	A2250
<b>Product Name:</b>	Dextran Sulfate 500 Sodium Salt <i>BioChemica</i>
<b>Specifications:</b>	average mol wt: approx. 500,000 g/mol Solubility (10 %; H <sub>2</sub> O): clear, colorless to yellow pH (1 %; H <sub>2</sub> O): 5.0 - 7.5 Sulfated ash: 40 - 50 % Loss on drying: max. 10 % Chloride: max. 0.1 % Fe: max. 0.0002 %
<b>WGK:</b>	1
<b>Storage:</b>	2 - 8°C
<b>CAS:</b>	9011-18-1
<b>CS:</b>	39139000

**AppliChem GmbH**