



NECK DATA: USE R3/33 SPECIFICATION,  
1 1/4 TURNS OF FULL THREAD TO BE  
DEPRESSED 0.2 AT MOULD JOINT.

ALPHA NUMERIC CODE EMBOSSED  
ON BASE 0.6 DEEP AS SHOWN.



SERIES LETTER, '30' AND BEATSON  
CLARK LOGO EMBOSSED 0.4 DEEP



VIEW ON BASE SHOWING  
EMBOSSED MARKING STANDING  
AREA EMBOSSED WITH BAR STIPPLE.

SCALE 1:1 (A2)		DESIGNER CH	30ml POWDER	
DATE 08/11/99	CHECKED DA			
NECK FINISH R3/33	REV A) NECK SPECIFICATION MODIFIED TO SUIT BS1918		DA 14/08/89.	NECK REF. No.
APPROX. WEIGHT 48 g	REV B) BOTTOM MARKING MODIFIED, DOT CODE DETAIL ADDED.		(DNA) DA 06/07/12.	<b>2A7/A1</b>
FILL POINT 13 down	REV C) KNURLING STIPPLE CHANGED TO BAR STIPPLE, ALPHA NUMERIC CODE ADDED. BASE PLATE DETAIL REPOSITIONED.		RL 25/04/2014.	JOB No.
CAPACITY TO FILL POINT 37 ml	<b>2985</b>			
CAPACITY TOLERANCE (BULK) ± 2.0 ml	DRAWING No.			
CAPACITY BRIMFULL 43 ml	<b>B3139/C</b>			
Note: Filling point measured from top centrally inside neck unless otherwise specified.	Untoleranced dimensions are for mould making guidance only and should be regarded as approximate.	This drawing and the design of the article in the drawing are each the property of . Neither the drawing nor the design shall be copied without the previous written authority of .		
DIMENSIONS IN MILLIMETRES				