


PDS No. 68261x / 682625 / 682660	<b>PRODUCT DATA SHEET</b>				Page 1 of 1
Revision 01	Roller Bottle CELLMASTER, 1 XL, PS				
	Greiner Item-No. 68261x / 682625 / 682660				
Valid for Item-No.:	<b>682612</b> (sterile)	<b>682615</b> (sterile)	<b>682625</b> (sterile)	<b>682660</b> (sterile)	

1.	Description / Specification	
1.1	Description	Roller Bottle, long, smooth surface, graduation (scale 100-2000 ml), physical surface treatment, sterile, for adherent cells <i>682612, -060: standard screw cap</i> <i>682615, -025: filter screw cap</i>
1.2	Dimensions	See Customer Drawing
1.3	Volume	Max. volume: 4970 ml Recommended working Volume: max. 450 ml Growth area: 1700 cm <sup>2</sup>
1.4	Material / Resin	<u>Bottle</u> : PS (Polystyrene), free of heavy metal <u>Screw cap</u> : HDPE (High Density Polyethylene), free of heavy metal <u>Filter</u> : PET (Polyethylene Terephthalate), free of heavy metal
1.5	Colour	<u>Bottle</u> : clear; print of graduation: black; writing area: white <u>Standard and filter screw cap</u> : blue <u>Filter</u> : white
1.6	Sterilisation	SAL 10 <sup>-3</sup>
1.7	Quality Control	- <u>Raw Material-Control</u> : physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for adherent cell culture to be used by qualified personnel in a laboratory environment.
1.9	Other Information	- For single use only - Expiry date and Lot-No. printed on label of bottle

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic. <i>682615, -625: controlled gas exchange through filter cap</i>
2.2	Temperature range	For application: -20°C to +60 °C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: <a href="https://www.gbo.com/en_INT/know-how-services/download-center.html">https://www.gbo.com/en_INT/know-how-services/download-center.html</a>
2.6	Shelf life	For application: 4 years
2.7	Other Information	<u>Bottle and cap</u> : USP Class VI tested

3.	Packaging	682612, -015	682625	682660
3.1	Pieces / Bag	12	12	1
3.2	Pieces / Box	12	24	12
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)		
3.4	Other Information	Certificate of Quality to download		

4.	Other Information
4.1	Laboratory use only

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 00	Date 18 December 2020	Date 15 January 2021	Date 25 February 2021	
Date 16.09.2020	Name S. Kaelberer	Name Dr. T. Schreiber	Name A. Illig	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.