



# DECLARATION DE CONFORMITE DECLARATION OF CONFORMITY



**FABRICANT  
MANUFACTURER**

**CORNING GOSSELIN SAS**

123 Rue de CAESTRE - CS 40019 BORRE  
59529 HAZEBROUCK Cedex - FRANCE  
Ph: +33.03.28.41.93.03, Fax: +33.03.28.49.56.92

**31/03/20**

**Gosselin, fabricant, déclare sous sa responsabilité que le(s) produit(s) :**  
**Manufacturer Gosselin, declares under his responsibility that the product(s):**

Référence	Designation
BTP30C-002	Container 40mL PP RedScrewCap Sterile 1/Bag 150/CS
BTP35C-002	Container 60mL PP RedScrewCap Sterile 1/Bag 100/CS
BTP52C-002	Container 125mL PP RedScrewCap Sterile 1/Bag 50/CS
BTP53C-002	Container 180mL PP RedScrewCap Sterile 1/Bag 50/CS
F2700-02	Bottle Square HDPE 2.7 L Handle YellowCapWad Assembled 30/CS
F2700-06	BottleSquareHDPE 2.7L Handle YellowRedCapWadAssembled30CS
F2700-12	BottleSquareHDPE2.7LBrown Handle YellowCapWad Assembled 30CS
LC2000-03	Bottle Square HDPE 2L 58mmRedCapSeal Assembled Sterile 50/CS
LC2000-06	Bottle Square HDPE 2L 58mm WhiteCap Seal Assembled 50/CS
LC2000-26	BottleSquare HDPE 2L 24HUrines RedCapSeal Nonassembled 50/CS
LC2000-32	Bottle Square HDPE 2L Label 58mmRed CapSeal Assembled 50/CS
LC2500A-01	Bottle Square HDPE 2.5 L 77mmWhiteCap Wad Assembled 35/CS
LC2500A-02	Bottle Square HDPE 2.5 L 77mmWhiteCap Wad Non-assembled 35/CS
LC2500A-03	Bottle Square HDPE 2.5L Brown 77mmWhiteCapWad Assembled 35CS
LC2500A-04	Bottle Square HDPE 2.5L Handle WhiteCapWad Assembled 35/CS
LC2500A-07	Bottle Square HDPE 2.5L 77mmWhiteCapWad AssembledSterile35CS
LC2500A-08	BottleSquareHDPE 2.5L Brown Handle WhiteCapWad Assembled35CS
LC2500A-09	BottleSquareHDPE 2.5L WhiteCapWad Assembled Sterile 1Bag35CS
LC500-25	BottleSquareHDPE 500mL LabelADDIS RedCapSeal Assembled 175CS
TP30C-002	Container 40mL PP RedScrewCap Sterile 100/Bag 1000/CS
TP30C-003	Container 40mL PP WhiteScrewCap Sterile 100/Bag 1000/CS
TP30C-004	Container 40mL PP GreenScrewCap Sterile 100/Bag 1000/CS
TP30C-011	Container 40mL PP BlueScrewCap 100/ Bag 1000/CS
TP30C-012	Container 40mL PP RedScrewCap 100/Bag 1000/CS
TP30C-013	Container 40mL PP WhiteScrewCap 100 /Bag 1000/CS
TP30C-014	Container 40mL PP GreenScrewCap 100 /Bag 1000/CS
TP30C-021	Container 40mL BluePP Label BlueScrewCap Sterile 100B 1000CS
TP30C-032	Container 40mL PP Label RedScrewCap Sterile 1000/CS
TP35C-001	Container 60mL PP BlueScrewCap Sterile 70/Bag 700/CS
TP35C-002	Container 60mL PP RedScrewCap Sterile 70/Bag 700/CS
TP35C-004	Container 60mL PP GreenScrewCap Sterile 70/Bag 700/CS
TP35C-011	Container 60mL PP BlueScrewCap 70/Bag 700/CS
TP35C-012	Container 60mL PP RedScrewCap 70/Bag 700/CS
TP35C-013	Container 60mL PP WhiteScrewCap 70/ Bag 700/CS

**Est(sont) conforme(s) avec :**

Les annexes I et III de la Directive européenne 98/79/CE du 27/10/1998 relative aux dispositifs médicaux de diagnostic in vitro. Ces dispositifs n'entrent pas dans la liste A ni dans la liste B de l'annexe II. Transposition française par l'ordonnance n°2001-198 du 1er/03/2001 et décret n°2004-108 du 04/02/2004.

*Cette déclaration est basée sur la documentation technique en vigueur démontrant la conformité du produit aux exigences de la directive.*

**Conform(s) with :**

Annex I & III of EC Directive 98/79/EC "In vitro" diagnostics Medical Devices adopted the 27th/10/1998.

*This declaration is based on the technical documentation showing the conformity of the product to the requirements of the directive*

*Signature Electronique / Electronic Signature*

**Cédric MAES**  
**Quality Assurance Manager**





# DECLARATION DE CONFORMITE DECLARATION OF CONFORMITY



**FABRICANT  
MANUFACTURER**

**CORNING GOSSELIN SAS**

123 Rue de CAESTRE - CS 40019 BORRE  
59529 HAZEBROUCK Cedex - FRANCE  
Ph: +33.03.28.41.93.03, Fax: +33.03.28.49.56.92

**31/03/20**

**Gosselin, fabricant, déclare sous sa responsabilité que le(s) produit(s) :**  
**Manufacturer Gosselin, declares under his responsibility that the product(s):**

Référence	Designation
TP35C-014	Container 60mL PP GreenScrewCap 70/ Bag 700/CS
TP35C-021	Container 60mL BluePP Label BlueScrewCap Sterile 70/B 700/CS
TP35C-025	Container 60mL WhitePP Spoon RedScrewCap Sterile 70/B 700/CS
TP35C-032	Container 60mL PP Label RedScrewCap Sterile 70/Bag 700/CS
TP35C-042	Container 60mL PP RedScrewCap Sterile 1/Bag 500/CS
TP35C-052	Container 60mL PP RedScrewCap Sterile 30/Tray 600/CS
TP35C-062	Container 60mL PP WhiteLabel RedCap Sterile 30/Tray 600/CS
TP52C-001	Container 125mL PP BlueScrewCap Sterile 380/CS
TP52C-002	Container 125mL PP RedScrewCap Sterile 380/CS
TP52C-003	Container 125mL PP WhiteScrewCap Sterile 380/CS
TP52C-004	Container 125mL PP GreenScrewCap Sterile 380/CS
TP52C-011	Container 125mL PP BlueScrewCap 380 /CS
TP52C-012	Container 125mL PP RedScrewCap Assembled 380/CS
TP52C-013	Container 125mL PP WhiteScrewCap 380/CS
TP52C-014	Container 125mL PP GreenScrewCap 380/CS
TP52C-021	Container 125mL BluePP WhiteLabel BlueScrewCap Sterile 380CS
TP52C-022	Container 125mL WhitePP WhiteSpatula RedScrewCap 380/CS
TP52C-023	Container 125mL WhitePP Spatula Red ScrewCap Sterile 380/CS
TP52C-024	Container 125mL WhitePP RedScrewCap 380/CS
TP52C-025	Container 125mL WhitePP BlueScrewCap 380/CS
TP52C-032	Container 125mL PP WhiteLabel RedScrewCap Sterile 380/CS
TP53C-001	Container 180mL PP BlueScrewCap Sterile 264/CS
TP53C-002	Container 180mL PP RedScrewCap Sterile 264/CS
TP53C-003	Container 180mL PP WhiteScrewCap Sterile 264/CS
TP53C-004	Container 180mL PP GreenScrewCap Sterile 264/CS
TP53C-011	Container 180mL PP BlueScrewCap 264 /CS
TP53C-012	
TP53C-013	Container 180mL PP WhiteScrewCap 264/CS
TP53C-014	Container 180mL PP GreenScrewCap 264/CS
TP53C-021	Container 180mL BluePP WhiteLabel BlueScrewCap Sterile 264CS
TP53C-023	Container 180mL WhitePP RedScrewCap Sterile 264/CS
TP53C-032	Container 180mL PP WhiteLabel RedScrewCap Sterile 264/CS
TP53C-042	Container 180mL PP RedScrewCap Sterile 1/Bag 200/CS
TP53C-052	Container 180mL PP WhiteSpatula Red ScrewCap Sterile 264/CS

**Est(sont) conforme(s) avec :**

Les annexes I et III de la Directive européenne 98/79/CE du 27/10/1998 relative aux dispositifs médicaux de diagnostic in vitro. Ces dispositifs n'entrent pas dans la liste A ni dans la liste B de l'annexe II.  
Transposition française par l'ordonnance n°2001-198 du 1er/03/2001 et décret n°2004-108 du 04/02/2004.

*Cette déclaration est basée sur la documentation technique en vigueur démontrant la conformité du produit aux exigences de la directive.*

**Conform(s) with :**

Annex I & III of EC Directive 98/79/EC "In vitro" diagnostics Medical Devices adopted the 27th/10/1998.

*This declaration is based on the technical documentation showing the conformity of the product to the requirements of the directive*

**Signature Electronique / Electronic Signature**

**Cédric MAES**  
**Quality Assurance Manager**

