


TDS No 0072		TECHNICAL DATA SHEET				
Revision 001		Microtitre Plates, Coloured				
Product Codes		611F96BK		611F96WT		
Colour		Black		White		
Well Capacity (µl)		400		400		
Well Format		Flat Bottom		Flat Bottom		
		Hanging Well		Chimney Well		
1	Dimensions (mm)	Length (1-12)	Width (A-H)	Pitch Well A1	Pitch Well A12	
		127.76	85.48	14.38	14.38	
2	Regulatory	CE	96 well microtitre plates meet the Essential requirements of the Council Directive, 98/79/EC 27th October 1998 Class 1 Devices.			
		SBS	Coloured microtitre plates conform to the dimensions specified in the Society for Biomolecular Screening standard 2006			
3	Material composition All materials are FDA approved		Virgin polystyrene			
4	Sterilisation	Non Sterile	NS			
5	Autoclavability		Not Suitable			
6	Chemical Resistance		Dilute Acids, Alcohols, Bases, Oils & Minerals			
7	Shelf Life		5 years from date of manufacture			
8	Packaging Format		Product Code	Shelf Pack Qty	Case Qty	
			611F96BK	5	50	
			611F96WT	5	50	
9	Other information		All Materials are non-cytotoxic			
Data Sheet subject to change without notice.						
Prior Issue		Raised	Approved	Released	Information contained within this document is confidential and proprietary to ThermoFisher Scientific. This document must not be reproduced. All rights of design, and copyright are reserved.	
Revision: N/A		Date:	Date:	Date:		
Date: N/A		Name: Sian Hayman	Name:	Name:		