

2mag SINGLE POINT STIRRER FOR CELL CULTURES AND TISSUE CULTURES

Designed for chemical, biotechnical and medical applications of all kinds, developed for maximum requirements, maximum mechanical, thermal and chemical stability as well as uncompromising continuous operation.



2mag bioMIX 1

2mag
magnetic (e)motion

Schragenhofstr. 35
80992 München

Tel. +49 (89) 14 33 42 52
Fax +49 (89) 14 33 43 69

magnetic **e**motion

Type	bioMIX 1
Stirring points	1
Stirring volume	5 - 5.000 ml
Speed range	5 - 250 rpm
Stirring power (max.)	6 W
Power setting	10-steps (10 – 100%)
Material housing	stainless-steel
Permitted operation conditions	-10 up to +50 °C (at 80% humidity)
Permitted storage conditions	-40 °C up to +70 °C, 10 – 95 %, 500 - 1060 hPa
Electrical data	100–240 volts / 50-60 Hz / 1.5 A
Measurements (BxTxH)	188 x 238 x 38 mm
Protection category	IP 64
Weight	4,0 kg
Order No.	85001

Individual units and dimensions on request

3 years warranty on material and manufacturing. We reserve the right to make technical alterations. We do not assume liability for errors in printing or typing.