

DATASHEET
PART NUMBER 760207

ClearLine ThermalSeal sealing films

Application :


PCR, long-term storage

Advantage :


Easy to remove

Compatible with :

Skirted and semi-skirted


Features

Packaging	Films
End-tabs	Features two end-tabs with perforations
Collar	Without collar
Material	Polypropylene
Adhesive material	Acrylic
Piercable	No
Chemically resistant	Yes
Breathable	No
Opacity	Transparent
Colour	-
Sterile	Yes - DNase RNase free
Functional temperature range	-40°C to +120°C
Film/Foil thickness	50,8 µm
Adhesive thickness	24,1 µm
Dimensions (without end-tabs)	123,1 x 79,4


Units : 100 films