

ClearLine Ezcap 1x8 PCR sealing films

Application :

PCR

Advantage :


Contains 6 strips on a continuous liner

Compatible with :

8-Tube Strip and Strip wells - Skirted and semi-skirted

Features

Packaging	Films
End-tabs	Features two end-tabs with perforations
Collar	Without collar
Material	Polyester
Adhesive material	Acrylic
Piercable	No
Chemically resistant	No
Breathable	No
Opacity	Transparent
Colour	-
Sterile	No - DNase RNase free
Functional temperature range	-20°C to +120°C
Film/Foil thickness	50,8 µm
Adhesive thickness	41,9 µm
Dimensions (without end-tabs)	73,3 x 10,1



Units : 400 films