





## SCHEDA TECNICA PRODOTTO - TECHNICAL DATA SHEET

### VIALS FOR DROSOPHILA



#### Characteristics

- In carton-rack of 100 pcs or bulk



- In polystyrene
- 2 formats of tubes:  
Ø 25 x 95 mm  
Ø 28.5 x 95 mm

Cat.no	Description and Dimension	Packaging (pcs)
789008	Drosophila tube Ø 25 x 95 mm in PS	500 bulk
789001	Drosophila tube Ø 25 x 95 mm in PS	5 tray x 100
789009	Drosophila tube Ø 28.5 x 95 mm in PS	500 bulk
789002	Drosophila tube Ø 28.5 x 95 mm in PS	5 tray x 100