

Description

- Droso-Plugs® Drosophila Closures

Part number

- 789193
- Pieces by sales : 6760 caps



Description

- Droso-Plugs® Drosophila Closures
- Fits Narrow Vials
- Densest foam ever offered to fly labs
- Reusable - washable & autoclavable
- Rec. Autoclave Settings: 15min, 121° & 15psi (avoid dry cycle)
- Remove Droso-Plug from vial/bottle before autoclaving
- Height: 1.125in. (28mm)
- Diameter: 0.984in. (25mm)
- The Droso-Plugs® are the densest foam closures available to fly labs. The proprietary nominal density of 92.1kg/m³ provides an incredibly tight cell structure (3,000,000 micro-cells per square centimeter). They are reusable (washable, autoclavable) making them the easiest closure ever available: Insert/Remove/Gas. Results may vary after autoclaving. Test a small amount first.

