

Corning® Specialty Media and Supplements Product Portfolio

CORNING



Corning offers a broad portfolio of Specialty Media and Supplements to meet the need of researchers and those requiring specialty and custom options for scale-up requirements. The variety of specialized media covers those who require serum-free and low-protein products, and to more specialized products such as for insect cell culture and monoclonal antibody production.

Corning specialty media and supplements are developed and commercialized for the optimum cultivation of specific cell types. Our product line offers several defined, serum-free media that offer numerous advantages over the use of serum, including reduced foaming, the elimination of exogenous agents, easier downstream processing.



Specialty Media


| Cat. No. | Description | Unit Size | Qty/Pk |
|--------------|---|-----------|--------|
| 355056 | Corning Hepatocyte Culture media kit | 500 mL | 1 |
| 25-072-CV | Lymphocyte Separation Medium (LSM)* | 500 mL | 1 |
| 40-101-CV | SF (Serum-free) medium, with L-glutamine and 1 g/L BSA | 500 mL | 6 |
| 40-102-CV | SF medium, with L-glutamine and 1 g/L BSA, without phenol red | 500 mL | 6 |
| 40-200-CV | PF (Protein-free) medium, without L-glutamine | 500 mL | 6 |
| 40-215-CVR | Corning hybrigro™ SF medium | 500 mL | 6 |
| 40-300-CVR | Corning transfectagro™ reduced-serum medium, with HEPES, Corning glutagro™, and sodium pyruvate | 500 mL | 6 |
| 46-060-CM | TSB (Tryptic Soy Broth), animal-free, liquid* | 1L | 1 |
| 355054 | Corning Endothelial Cell Culture medium kit | 500 mL | 1 |
| 40-05-100-1A | Corning NutriStem® hPSC SF medium* | 500 mL | 1 |
| 355058 | Corning Intestinal Differentiation media pack, kit | – | 1 |

*Other varieties are also available.



Insect Cell Media

| Cat. No. | Description | Unit Size | Qty/Pk |
|-----------|--|-----------|--------|
| 13-402-CV | Corning insectagro® DS2 serum-free/protein-free medium, 1X without L-glutamine | 500 mL | 6 |
| 13-410-CV | Corning insectagro Sf9 serum-free/protein-free medium, 1X [+]-L-glutamine | 500 mL | 6 |



Corning's islet solutions were designed and tailored specifically for diabetes research applications. A complete specialty line for pancreatic islet isolation and maintenance is also available.

For more information, visit www.corning.com/lifesciences.

Corning also sells a whole line of antibiotics and specialty supplements.



| Cat. No. | Description | Unit Size | Qty/Pk |
|-----------|------------------------------------|-----------|--------|
| 61-242-RG | Tetracycline hydrochloride, powder | 5g | 1 |
| 30-005-CR | Gentamicin sulfate, liquid | 10 mL | 1 |
| 61-098-RF | Gentamicin sulfate, powder* | 100 mg | 1 |
| 30-240-CR | Hygromycin B solution | 10 mL | 1 |
| 30-234-CR | G418 sulfate, liquid* | 50 mg/mL | 1 |
| 61-234-RF | G418 sulfate, powder* | 1 g | 1 |
| 30-006-CF | Kanamycin sulfate, liquid | 50 mL | 1 |
| 61-176-RG | Kanamycin sulfate, powder* | 5 g | 1 |
| 61-241-RG | Neomycin sulfate, powder | 5 g | 1 |
| 61-238-RH | Ampicillin, sodium salt, powder* | 10 g | 1 |

* Other varieties are also available.

Corning® Media Growth Supplements

| Cat. No. | Description | Unit Size | Qty/Pk |
|----------|--|-----------|--------|
| 354123 | Bovine pituitary extract (BPE) | 15 mg | 1 |
| 356006 | Endothelial cell growth supplement (ECGS) | 100 mg | 1 |
| 354351 | ITS premix universal culture supplement | 5 mL | 1 |
| 354352 | ITS+ premix universal culture supplement | 20 mL | 1 |
| 355006 | MITO+ serum extender | 5 mL | 1 |
| 355100 | Corning Nu-Serum growth medium supplement | 100 mL | 1 |
| 355104 | Corning Nu-Serum IV growth medium supplement | 100 mL | 1 |
| 354115 | T-cell culture supplement with ConA (IL-2 culture supplement), rat | 100 mL | 1 |

For more specific information on claims, visit the Certificates page at www.corning.com/lifesciences.

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

For additional product or technical information, visit www.corning.com/lifesciences or call 800.492.1110. Outside the United States, call +1.978.442.2200 or contact your local Corning sales office.

CORNING

**Corning Incorporated
Life Sciences**

836 North St.
Building 300, Suite 3401
Tewksbury, MA 01876
t 800.492.1110
t 978.442.2200
f 978.442.2476

www.corning.com/lifesciences

**ASIA/PACIFIC
Australia/New Zealand**

t 61 427286832
China
t 86 21 3338 4338
f 86 21 3338 4300

India
t 91 124 4604000
f 91 124 4604099

Japan
t 81 3-3586 1996
f 81 3-3586 1291

Korea
t 82 2-796-9500
f 82 2-796-9300

Singapore
t 65 6572-9740
f 65 6861-2913

Taiwan
t 886 2-2716-0338
f 886 2-2516-7500

**EUROPE
CSEurope@corning.com**

France
t 0800 916 882
f 0800 918 636

Germany
t 0800 101 1153
f 0800 101 2427

The Netherlands
t 020 655 79 28
f 020 659 76 73

United Kingdom
t 0800 376 8660
f 0800 279 1117

**All Other European
Countries**
t +31 (0) 206 59 60 51
f +31 (0) 206 59 76 73

**LATIN AMERICA
grupoLA@corning.com**

Brasil
t 55 (11) 3089-7400

Mexico
t (52-81) 8158-8400