

Certificate of Analysis

Description: Tris Base Buffer

(MW121.14)

CAS # 77-86-1

Storage: 15-30°C.

Catalog #: 61-233-RR

Lot #: 61233094

Expiration Date: 2017/11

Country of Origin/Final Processing: USA

Notes:

This product was tested in accordance with currently approved Mediatech, Inc. specifications and procedures. Testing procedures are maintained in compliance with the current versions of the USP and/or EP, where applicable.

For Research Use Only. Not for use in Diagnostic or Therapeutic procedures. Utilization of this product apart from the labeled intended use may be a violation of local and/or Federal Law.

TEST PARAMETERS	SPECIFICATION	RESULT
Appearance	White crystals/crystalline powder	Passes
DNase	None Detected	None Detected
RNase	None Detected	None Detected
Protease	None Detected	None Detected

For detailed information on Mediatech, Inc. testing procedures, please contact our Technical Services department at (800) 235-5476 or via email at MTTechServ@corning.com.

Written By/Date:

J. Blanco 01 Oct 2014
Quality Systems Department

Verified By/ Date:

[Signature] 02 Oct 2014