



Performance You Can Trust Quality, Accuracy, Cleanliness



Reputation for quality

Precision-molded pipette tips from Rainin have an established reputation for accuracy, quality and cleanliness within labs worldwide. Rainin FinePoint™ tip design and stringent quality testing provide highly accurate pipette tips.




Certified clean tips

Tips are manufactured in Class 100,000 cleanroom conditions. Virgin polypropylene material is fed from silos directly into molding machines, preventing human contact. Each lot of tips is rigorously tested and certified to be free of detectable biological contaminants



Hinged-lid rack

Open the racks with one hand – simply press the latch and lift the lid. Racks stack to save space. Colored tip decks make the volume readily apparent and allow full or partial loading of multichannel pipettes.



Rainin tip filters

Unlike other tip filters, Rainin filters contain no carboxymethyl cellulose (CMC), which can disrupt PCR analyses. They consist of polyethylene granules that have been sintered into a three-dimensional structure containing numerous pores and convoluted channels. These filters arrest aerosol particles from the airstream as they impact against the interior of the filter.



Research Grade Pipette Tips Exceptional Quality and Value

Rainin Research Grade pipette tips are designed to achieve a comfortable, low-force fit on most major pipette brands, including Rainin™, Gilson®, BioHit®, Eppendorf® and ThermoFisher®.

Visual tip inspection and real-time quality testing during production ensure the absence of defects that could potentially compromise tip performance. Class 100,000 cleanroom manufacturing prevents tip contamination during molding and packaging.

Rainin Research Grade Pipette Tips

For easy, accurate pipetting, Research Grade tips fit all major brand pipettes

- **High quality and purity**

Research Grade tips are certified to be free of biological contamination and come with a certificate of compliance.

- **Filter option to prevent cross contamination**

Rainin filters are made of pure sintered polyethylene to prevent aerosols from contaminating the pipette.

- **Extended length available for tall, narrow vessels**

10 µL extended-length pipette tips are available in bulk and rack formats: non-sterile, sterile and filtered.



Ordering information

MT-Order No.	Description	Max. Capacity
Racked Tips		
30374654	Tips HRC UNV 10µl 960C/10	960 tips in 10 racks
30374657	Tips HRC UNV 10µl S 960C/10	960 tips in 10 racks, Sterilized
30374660	Tips HRC UNV 10µl F 960C/10	960 tips in 10 racks, Sterilized, Filter
30537525	Tips HRC UNV 10µl X 960C/10	960 tips in 10 racks, Extended
30537526	Tips HRC UNV 10µl SX 960C/10	960 tips in 10 racks, Extended, Sterilized
30537527	Tips HRC UNV 10µl XF 960C/10	960 tips in 10 racks, Extended, Sterilized, Filter
30374661	Tips HRC UNV 20µl F 960C/10	960 tips in 10 racks, Sterilized, Filter
30537521	Tips HRC UNV 100µl F 960C/10	960 tips in 10 racks, Sterilized, Filter
30537522	Tips HRC UNV 200µl F 960C/10	960 tips in 10 racks, Sterilized, Filter
30374655	Tips HRC UNV 250µl 960C/10	960 tips in 10 racks
30374658	Tips HRC UNV 250µl S 960C/10	960 tips in 10 racks, Sterilized
30374656	Tips HRC UNV 1000µl 768C/8	768 tips in 8 racks
30374659	Tips HRC UNV 1000µl S 768C/8	768 tips in 8 racks, Sterilized
30537523	Tips HRC UNV 1000µl F 768C/8	768 tips in 8 racks, Sterilized, Filter
RSC Tips		
30537471	Tips RSC UNV 10µl 960C/10	960 tips in 10 racks
30537474	Tips RSC UNV 10µl S 960C/10	960 tips in 10 racks, Sterilized
30537472	Tips RSC UNV 250µl 960C/10	960 tips in 10 racks
30537475	Tips RSC UNV 250µl S 960C/10	960 tips in 10 racks, Sterilized
30537473	Tips RSC UNV 1000µl 768C/8	768 tips in 8 racks
30537476	Tips RSC UNV 1000µl S 768C/8	768 tips in 8 racks, Sterilized
Bulk Tips		
30537470	Tips RC UNV 10µl 1000C	1000 tips in a small bag
30537524	Tips RC UNV 10µl X 1000C	1000 tips in a small bag, Extended
30537469	Tips RC UNV 250µl 1000C	1000 tips in a small bag
30537468	Tips RC UNV 1000µl 1000C	1000 tips in a small bag

Rainin, Pipetting 360° and FinePoint are trademarks of Mettler-Toledo Rainin, LLC. Gilson, BioHit, Eppendorf and ThermoFisher are trademarks of their respective owners.

Mettler-Toledo International Trading (Shanghai) Co., Ltd
 589 Guiping Road
 Shanghai, 200233, P.R. China
 Tel 86-21-64850435
 Fax 86-21-64853351

Mettler-Toledo (ChangZhou) Measurement Technology Ltd.
 111 West Taihu road, Xinbei District
 Changzhou, Jiangsu, 213125 P.R. China
 Tel 86-519-86642040
 Fax 86-519-86641991 17701488 Rev. A

www.mt.com/rainin

For more information