

## TECHNICAL DATASHEET

## PART NUMBER

**713249**




The ISO8655 norm guarantees that the tips used do not alter the technical characteristics of the pipettes, as given by the manufacturer.



### ClearLine Tips

- ✓ In clear polypropylene
- ✓ Store in a clean and dry place, at room temperature
- ✓ Each batch number is certified : RNase, DNase free, Human DNA free, PCR inhibitor free and pyrogen free
- ✓ Traceability : 5 digit batch number, part number, volumetric tips, volumetric tips and their expiration date.

|               |                                     |
|---------------|-------------------------------------|
| Part number   | 713249                              |
| Packaging     | ClearLine plus rack with hinged lid |
| Volume        | 200 µl                              |
| Graduation    | 10, 50, 100                         |
| Length        | 51,44 mm                            |
| Color         | Natural                             |
| Filter        | No                                  |
| Low Retention | No                                  |
| Sterility     | Yes                                 |
| Autoclavable  | Yes                                 |



### Sales units :

960 tips