


TECHNICAL DATASHEET
PART NUMBER
713254



The ISO8655 norm guarantees that the tips used do not alter the technical characteristics of the pipettes, as given by the manufacturer.



ClearLine Tips

- ✓ In clear polypropylene
- ✓ Store in a clean and dry place, at room temperature
- ✓ Each batch number is certified : RNase, DNase free, Human DNA free, PCR inhibitor free and pyrogen free
- ✓ Traceability : 5 digit batch number, part number, volumetric tips, volumetric tips and their expiration date.

| | |
|----------------------|-------------------------------------|
| Part number | 713254 |
| Packaging | ClearLine plus rack with hinged lid |
| Volume | 1250 µl |
| Graduation | 200, 500, 1000 |
| Length | 102 mm |
| Color | Natural |
| Filter | No |
| Low Retention | No |
| Sterility | No |
| Autoclavable | Yes |


Sales units :

768 tips