

External dimensions

Height (mm)	437
Width (mm)	588
Depth(mm)	288

Internal dimensions

Height (mm)	245
Width (mm)	460
Depth (mm)	180

Volumes

Gross volume (litres)	20
-----------------------	----

Ship dimensions

Weight net (kg)	18
Height (mm) Cardboard	470
Width (mm) Cardboard	610
Depth (mm) Cardboard	310
Volume (m3) Cardboard	0.089
Weight (kg) Cardboard	19

Height (mm) Wood	-
Width (mm) Wood	-
Depth (mm) Wood	-
Volume (m3) Wood	-
Weight (kg) Wood	-

Accessories

12 Ice Packs PCM -32°C of 0.6L Metal Box and lid

Option

Testo Logger 174T/ starter set or LogTag TRID30-7R or LogTag TRED30-7R
--

Passive Blood transport box MT8-ET-32

Technical details

Colour	blue according to RAL5015
CFC/HCFC free	yes
External / internal material	Polyethylene
Insulation material	50-60 mm Polyurethane foam / Foaming agent; Cyclopentane
Weight fully loaded (kg)	18.91
Required Ice Packs	-
Required Eutetics	12 x 0.6L at ambient +22°C 12 x 0.6L at ambient +43°C

Performance

Hold over time at -10°C	> 24h
Hold over time at +22°C	> 24h
Hold over time at +32°C	-
Hold over time at 43°C	> 24h
Cold life without lid openings +2°C to +10°C at ambient 32°C	-
Cold life without lid openings +2°C to +10°C at ambient 43°C	-
Cold life without lid openings -3°C to +10°C	-
Cold life at +43°C	-
Warm life at -20°C	-
Cool life at +43°C	-
Blood storage capacity (litres)	6 bags of 300ml
Climate Zone	-

Certification

WHO test procedure	-
WHO specifications	-
WHO test report	-
PQS Code	-
PIS Code	-
ADR/RID MARITIM IATA / ICAO	UN 4H2 / Y19 / S / "YYYY" B / 921-110058
Safety standards	NA
EMC standards	NA
Medical Device Regulation	(EU) 2017/745: Class IIa