

TECHNICAL DATA SHEET

Part Number	71696-2.5L-HON
Description	Sodium hypochlorite solution

Technical specifications :

- ✓ sodium hypochlorite solution, purum, ~10% (rt) 2,5l
- ✓ CAS N°7681-52-9
- ✓ Chemical purity $\geq 10\%$
- ✓ Molecular weight 74.44 g/mol
- ✓ Molecular formula ClNaO
- ✓ Density 1.206 g/cm³ (25 °C)
- ✓ Flash point Not applicable
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Bottle
- ✓ Package size 2.5L

- ✓ GHS pictograms GHS09 ; GHS05
- ✓ H Statements H290 ; H314 ; H410
- ✓ P Statements P284 ; P280 ; P260

Sales unit :

2.5L Sodium hypochlorite solution
