



Honeywell Specialty Chemicals Seelze GmbH
Wunstorferstrasse 40
Seelze, 30926
Telefon: +49 5137 999-0
lab.honeywell.com

Brand: Honeywell
Product: 71890
CAS Numbers (All) : 7775-27-1
Molar Mass : 238.1 g/mol

Sodium persulfate purum p.a., 99.0% (RT)

Parameter	Specification	Units
appearance	complying	
assay	min. 99.0	%
calcium (Ca)	max. 50	ppm
copper (Cu)	max. 50	ppm
cadmium (Cd)	max. 50	ppm
cobalt (Co)	max. 50	ppm
iron (Fe)	max. 50	ppm
potassium (K)	max. 100	ppm
nickel (Ni)	max. 50	ppm
lead (Pb)	max. 50	ppm
zinc (Zn)	max. 50	ppm
chloride (Cl)	max. 50	ppm