


PDS No. 73118x	PRODUCT DATA SHEET			Page 1 of 1
Revision 03	Disposable Inoculation Needle, PS			 greiner bio-one
	Greiner Item-No. 73118x			
Valid for Item-No.:	731180 (sterile)	731181 (sterile)	731185 (sterile)	

1.	Description / Specification	
1.1	Description	PS disposable inoculation needle
1.2	Dimensions	See Customer Drawing
1.3	Volume	-
1.4	Material / Resin	PS (Polystyrene), free of heavy metal
1.5	Colour	Yellow
1.6	Sterilization	SAL 10 ⁻³
1.7	Quality Control	- <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	Hydrophobic
2.2	Temperature range	-20°C to +60°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	5 years after month of production
2.7	Other Information	-

3.	Packaging	731180	731181	731185
3.1	Pieces / Bag	50	1	10
3.2	Pieces / Box	2000	600	3000
3.3	Lot-No.	F YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	F YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	F YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	-	-	

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 02	Date 21 April 2015	Date 22 April 2015	Date 22 April 2015	
Date 20.11.2006	Name S. Kaelberer	Name Dr. L. Breth	Name A. Schulz	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.