

ClearLine MiniStrip sealing films

Application :


Short-Term Storage, Incubation, ELISA

Advantage :


Contains 4 strips on a continuous liner

Compatible with :

8-Tube Strip and Strip wells - Skirted and semi-skirted


Features

Packaging	Films
End-tabs	Features two end-tabs with perforations
Collar	Without collar
Material	Polypropylene
Adhesive material	Acrylic
Piercable	No
Chemically resistant	No
Breathable	No
Opacity	Transparent
Colour	-
Sterile	No
Functional temperature range	-40°C to +120°C
Film/Foil thickness	50,8 µm
Adhesive thickness	26,7 µm
Dimensions (without end-tabs)	82,6 x 19,1


Units : 200 films