

### ClearLine MiniStrip sealing films

**Application :**


Short-Term Storage, Incubation, ELISA

**Advantage :**


Contains 4 strips on a continuous liner

**Compatible with :**

8-Tube Strip and Strip wells - Skirted and semi-skirted


**Features**

|                                      |                                         |
|--------------------------------------|-----------------------------------------|
| <b>Packaging</b>                     | Films                                   |
| <b>End-tabs</b>                      | Features two end-tabs with perforations |
| <b>Collar</b>                        | Without collar                          |
| <b>Material</b>                      | Polypropylene                           |
| <b>Adhesive material</b>             | Acrylic                                 |
| <b>Piercable</b>                     | No                                      |
| <b>Chemically resistant</b>          | No                                      |
| <b>Breathable</b>                    | No                                      |
| <b>Opacity</b>                       | Transparent                             |
| <b>Colour</b>                        | -                                       |
| <b>Sterile</b>                       | Yes                                     |
| <b>Functional temperature range</b>  | -40°C to +120°C                         |
| <b>Film/Foil thickness</b>           | 50,8 µm                                 |
| <b>Adhesive thickness</b>            | 26,7 µm                                 |
| <b>Dimensions (without end-tabs)</b> | 82,6 x 19,1                             |


**Units : 200 films**