

ClearLine ThinSeal sealing films

Application :


Short-Term Storage, Incubation, ELISA

Advantage :


Marginally pierceable with single-channel pipettes for sample recovery

Compatible with :

Skirted and semi-skirted


Features

Packaging	Films
End-tabs	Features two end-tabs with perforations
Collar	Without collar
Material	Polyester
Adhesive material	Acrylic
Pierceable	Yes
Chemically resistant	No
Breathable	No
Opacity	Transparent
Colour	-
Sterile	Yes
Functional temperature range	-40°C to +120°C
Film/Foil thickness	25,4 µm
Adhesive thickness	24,1 µm
Dimensions (without end-tabs)	121,2 x 79,4


Units : 100 films