

ClearLine ThermalSeal RT2 sealing films

Application :


qPCR

Advantage :


Optically transparent sealing film with a very strong adhesive

Compatible with :

Two component plate -Skirted and semi-skirted


Features

Packaging	Films
End-tabs	Features two end-tabs with perforations
Collar	Without collar
Material	Polyester
Adhesive material	Acrylic
Piercable	No
Chemically resistant	No
Breathable	No
Opacity	Optically Clear
Colour	-
Sterile	No - DNase RNase free
Functional temperature range	-40°C to +120°C
Film/Foil thickness	50,8 µm
Adhesive thickness	29,2 µm
Dimensions (without end-tabs)	121,9 x 79,4


Units : 100 films